



Thermal Properties Analyzer

Liquids, Solids, Paste & Powders

ThermTest Inc
2 Garland Court, Suite 006
Fredericton, NB E3B 6C2
Canada

Phone: 506-458-5350
Fax: 866-274-5269
E-mail: info@ThermTest.com

www.ThermTest.com



Includes:

- KD2 Controller
- KD2 Needle sensor
- Conductive paste
- Carrying Case

Measurement Speed: 1.5 minutes

Accuracy: Thermal Conductivity 5%

- 0.02 to 2 W/mK
- -20° to 60° C

IEEE Standard 442-1981 and ASTM Standard D 5334-00.